



U.O.C. Neonatologia-TIN-Nido
Responsabile: Dott.^{ssa} Isabella Mondello

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9. M. Clementi. Effect of a chlorinated disinfectant (AMUCHINA, AMUCHINA S.p.A. Genova, Italy) on hepatitis C Virus (HCV) in vitro: Analysis of HCV binding to the cell surface receptors and analysis of viral replicators. *Acts Toxicol. Ther.* 18.1.1997

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11. Acute oral toxicity study in rats treated with the test article AMUCHINA Electrolytic Chloroxidizer, RBM, Torino, 1991.

12. Acute intravenous toxicity study in rats treated with the test article AMUCHINA Electrolytic Chloroxidizer RBM. Torino 1991

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Report No: ASH17-033560-03

Date: Dec 15 2017

Test time: The test was conducted at 24 hours post cycle completion.

Test Item	Test Group	Number of recovered microorganism (CFU/sample)	Efficacy(%)	Average Efficacy(%)	Log Reduction	Judge according to client's requirement: Log Reduction>6
Feeding bottle	1	<1	>99.9999	>99.9999	>6.1461	Pass
	2	<1	>99.9999			
Feet	1	<1	>99.9999	>99.9999	>6.1461	Pass
	2	<1	>99.9999			
Screw	1	<1	>99.9999	>99.9999	>6.1461	Pass
	2	<1	>99.9999			
Cover	1	<1	>99.9999	>99.9999	>6.1461	Pass
	2	<1	>99.9999			
Pump Valve	1	<1	>99.9999	>99.9999	>6.1461	Pass
	2	<1	>99.9999			
Pump Body	1	<1	>99.9999	>99.9999	>6.1461	Pass
	2	<1	>99.9999			

Remark:

Salmonella enterica subsp. Enterica ATCC 14028 was historically known as Salmonella typhimurium.

Test temperature: 24.7°C

Test humidity: 54.0%


SAMPLE DESCRIPTION: electric steam steriliser

*** End ***

SGS-CSTC Standards Technical Services (Shanghai) Co.,Ltd.

Page 9 of 9

Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com.cn/Terms-and-Conditions.aspx>, and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at [http://www.sgs.com.cn/Terms-and-Conditions.aspx](http://www.sgs.com.cn/Terms-and-Conditions/Terms-and-Conditions.aspx). Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full without the written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested.

Building No.889, Yishan Road, Shanghai, China 200233

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Report No: ASH17-033560-03


Date: Dec 15 2017

Test time: The test was conducted at 24 hours post cycle completion.

Test Item	Test Group	Number of recovered microorganism (CFU/sample)	Efficacy(%)	Average Efficacy(%)	Log Reduction	Judge according to client's requirement: Log Reduction>6
Feeding bottle	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Feet	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Screw	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Cover	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Pump Valve	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Pump Body	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			

Test microorganism: *Listeria monocytogenes* ATCC 7644
 Concentration of inoculated microorganism: 1.4×10^7 CFU/mL
 The initial bioburden: 1.4×10^6 CFU/sample
 Test time: The test was conducted at 0 hour post cycle completion.

Test Item	Test Group	Number of recovered microorganism (CFU/sample)	Efficacy(%)	Average Efficacy(%)	Log Reduction	Judge according to client's requirement: Log Reduction>6
Feeding bottle	1	<1	>99.9999	>99.9999	>6.1461	Pass
	2	<1	>99.9999			
Feet	1	<1	>99.9999	>99.9999	>6.1461	Pass
	2	<1	>99.9999			
Screw	1	<1	>99.9999	>99.9999	>6.1461	Pass
	2	<1	>99.9999			
Cover	1	<1	>99.9999	>99.9999	>6.1461	Pass
	2	<1	>99.9999			
Pump Valve	1	<1	>99.9999	>99.9999	>6.1461	Pass
	2	<1	>99.9999			
Pump Body	1	<1	>99.9999	>99.9999	>6.1461	Pass
	2	<1	>99.9999			

Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/terms-and-conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/terms-and-conditions/terms-e-documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that e-format content herein reflects the Company's findings at the time of its intervention only and with the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested.

Test Report

Report No: ASH17-033560-03


Date: Dec 15 2017

Test time: The test was conducted at 24 hours post cycle completion.

Test Item	Test Group	Number of recovered microorganism (CFU/sample)	Efficacy(%)	Average Efficacy(%)	Log Reduction	Judge according to client's requirement: Log Reduction>6
Feeding bottle	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Feet	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Screw	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Cover	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Pump Valve	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Pump Body	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			

Test microorganism: Cronobacter sakazakii ATCC 29004
 Concentration of inoculated microorganism: 1.1×10^7 CFU/mL
 The initial bioburden: 1.1×10^6 CFU/sample
 Test time: The test was conducted at 0 hour post cycle completion.

Test Item	Test Group	Number of recovered microorganism (CFU/sample)	Efficacy(%)	Average Efficacy(%)	Log Reduction	Judge according to client's requirement: Log Reduction>6
Feeding bottle	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Feet	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Screw	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Cover	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Pump Valve	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Pump Body	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			

Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/terms-and-conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/terms-and-conditions/terms-e-document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and with the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this Document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested.

Test Report

Report No: ASH17-033560-03

Date: Dec 15 2017

Test time: The test was conducted at 24 hours post cycle completion.

Test Item	Test Group	Number of recovered microorganism (CFU/sample)	Efficacy(%)	Average Efficacy(%)	Log Reduction	Judge according to client's requirement: Log Reduction>6
Feeding bottle	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Feet	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Screw	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Cover	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Pump Valve	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Pump Body	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			

Test microorganism: Streptococcus agalactiae ATCC 12386

Concentration of inoculated microorganism: 1.1×10^7 CFU/mL


The initial bioburden: 1.1×10^6 CFU/sample

Test time: The test was conducted at 0 hour post cycle completion.

Test Item	Test Group	Number of recovered microorganism (CFU/sample)	Efficacy(%)	Average Efficacy(%)	Log Reduction	Judge according to client's requirement: Log Reduction>6
Feeding bottle	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Feet	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Screw	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Cover	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Pump Valve	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			
Pump Body	1	<1	>99.9999	>99.9999	>6.0414	Pass
	2	<1	>99.9999			

SGS-CSTC Standards Technical Services (Shanghai) Co.,Ltd.

Page 6 of 9

Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service posted overleaf, available on request or accessible at <http://www.sgs.com/en/terms-and-conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/terms-and-conditions/terms-e-documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction clauses defined therein. Any holder of this document is advised that information contained herein reflects the Company's findings at the time of its intervention only and with the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not create any liability to a third party from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full without the written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested.

3F Building, No. 889, Yishan Road, Shanghai, China 200233

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Report No: ASH17-033560-03

Date: Dec 15 2017

Test time: The test was conducted at 24 hours post cycle completion.

Test Item	Test Group	Number of recovered microorganism (CFU/sample)	Efficacy(%)	Average Efficacy(%)	Log Reduction	Judge according to client's requirement: Log Reduction>6
Feeding bottle	1	<1	>99.9999	>99.9999	>6.2041	Pass
	2	<1	>99.9999			
Feet	1	<1	>99.9999	>99.9999	>6.2041	Pass
	2	<1	>99.9999			
Screw	1	<1	>99.9999	>99.9999	>6.2041	Pass
	2	<1	>99.9999			
Cover	1	<1	>99.9999	>99.9999	>6.2041	Pass
	2	<1	>99.9999			
Pump Valve	1	<1	>99.9999	>99.9999	>6.2041	Pass
	2	<1	>99.9999			
Pump Body	1	<1	>99.9999	>99.9999	>6.2041	Pass
	2	<1	>99.9999			

Test microorganism: Pseudomonas aeruginosa ATCC 9027

Concentration of inoculated microorganism: 1.2×10^7 CFU/mL


The initial bioburden: 1.2×10^6 CFU/sample

Test time: The test was conducted at 0 hour post cycle completion.

Test Item	Test Group	Number of recovered microorganism (CFU/sample)	Efficacy(%)	Average Efficacy(%)	Log Reduction	Judge according to client's requirement: Log Reduction>6
Feeding bottle	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Feet	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Screw	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Cover	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Pump Valve	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Pump Body	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			

SGS-CSTC Standards Technical Services (Shanghai) Co.,Ltd.

Page 5 of 9

Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com.cn/Terms-and-Conditions.aspx> and, for electronic format documents subject to Terms and Conditions for Electronic Documents printed at <http://www.sgs.com.cn/Terms-and-Conditions/Terms-and-Conditions-E-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction clauses defined therein. Any holder of this document is advised that information contained herein reflects the Company's findings at the time of its intervention only and with the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested.

1st Building, No. 889, Yishan Road, Shanghai, China 200233

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Report No: ASH17-033560-03

Date: Dec 15 2017

Test time: The test was conducted at 24 hours post cycle completion.

Test Item	Test Group	Number of recovered microorganism (CFU/sample)	Efficacy(%)	Average Efficacy(%)	Log Reduction	Judge according to client's requirement: Log Reduction>6
Feeding bottle	1	<1	>99.9999	>99.9999	>6.3222	Pass
	2	<1	>99.9999			
Feet	1	<1	>99.9999	>99.9999	>6.3222	Pass
	2	<1	>99.9999			
Screw	1	<1	>99.9999	>99.9999	>6.3222	Pass
	2	<1	>99.9999			
Cover	1	<1	>99.9999	>99.9999	>6.3222	Pass
	2	<1	>99.9999			
Pump Valve	1	<1	>99.9999	>99.9999	>6.3222	Pass
	2	<1	>99.9999			
Pump Body	1	<1	>99.9999	>99.9999	>6.3222	Pass
	2	<1	>99.9999			

Test microorganism: Salmonella enterica subsp. Enterica ATCC 14028

Concentration of inoculated microorganism: 1.9×10^7 CFU/mL


The initial bioburden: 1.9×10^6 CFU/sample

Test time: The test was conducted at 0 hour post cycle completion.

Test Item	Test Group	Number of recovered microorganism (CFU/sample)	Efficacy(%)	Average Efficacy(%)	Log Reduction	Judge according to client's requirement: Log Reduction>6
Feeding bottle	1	<1	>99.9999	>99.9999	>6.2041	Pass
	2	<1	>99.9999			
Feet	1	<1	>99.9999	>99.9999	>6.2041	Pass
	2	<1	>99.9999			
Screw	1	<1	>99.9999	>99.9999	>6.2041	Pass
	2	<1	>99.9999			
Cover	1	<1	>99.9999	>99.9999	>6.2041	Pass
	2	<1	>99.9999			
Pump Valve	1	<1	>99.9999	>99.9999	>6.2041	Pass
	2	<1	>99.9999			
Pump Body	1	<1	>99.9999	>99.9999	>6.2041	Pass
	2	<1	>99.9999			

SGS-CSTC Standards Technical Services (Shanghai) Co.,Ltd.

Page 4 of 9

Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/terms-and-conditions.aspx> and, for electronic format documents subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/terms-and-conditions/terms-e-document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction clauses defined therein. Any risk of this document is advised that information contained herein reflects the Company's findings at the time of its intervention only and with the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full without the written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested.

F Building, No 389, Yishan Road, Shanghai, China 200233

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Report No: ASH17-033560-03


Date: Dec 15 2017

Test time: The test was conducted at 24 hours post cycle completion.

Test Item	Test Group	Number of recovered microorganism (CFU/sample)	Efficacy(%)	Average Efficacy(%)	Log Reduction	Judge according to client's requirement: Log Reduction>6
Feeding bottle	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Feet	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Screw	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Cover	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Pump Valve	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Pump Body	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			

Test microorganism: Staphylococcus aureus ATCC 6538
 Concentration of inoculated microorganism: 2.1×10^7 CFU/mL
 The initial bioburden: 2.1×10^6 CFU/sample
 Test time: The test was conducted at 0 hour post cycle completion.

Test Item	Test Group	Number of recovered microorganism (CFU/sample)	Efficacy(%)	Average Efficacy(%)	Log Reduction	Judge according to client's requirement: Log Reduction>6
Feeding bottle	1	<1	>99.9999	>99.9999	>6.3222	Pass
	2	<1	>99.9999			
Feet	1	<1	>99.9999	>99.9999	>6.3222	Pass
	2	<1	>99.9999			
Screw	1	<1	>99.9999	>99.9999	>6.3222	Pass
	2	<1	>99.9999			
Cover	1	<1	>99.9999	>99.9999	>6.3222	Pass
	2	<1	>99.9999			
Pump Valve	1	<1	>99.9999	>99.9999	>6.3222	Pass
	2	<1	>99.9999			
Pump Body	1	<1	>99.9999	>99.9999	>6.3222	Pass
	2	<1	>99.9999			

Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com.cn/terms-sgs-conditions.aspx> and, for electronic documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com.cn/terms-and-conditions-terms-electronic-documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction clauses defined therein. Any holder of this document is advised that, if a formal or contented person reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not constitute a letter to a third party from assuming all their rights and obligations under the transaction documents. This document cannot be reproduced except in full without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested.

Test Report

Report No: ASH17-033560-03

Date: Dec 15 2017

TEST METHOD(S):

Method applied by client: In this study, the sterilization efficacy of Philips Avent 3-in-0 Electric Steam Steriliser (SCF284) will be evaluated at 0 hour as well as 24 hour poststerilization by determining the recovery of pre-loaded pathogen of 7 microorganism with a bioburden level of over 6 logs on different feeding appliance (feeding bottle, teat, pump valve, pump body, screw, cover). The ability to reduce microbial bioburden will be determined by comparing level of viable organisms recovered from the different treated feeding accessories to initial bioburden. The test will be conducted at both 0 hour and 24 hours post cycle completion.

TEST ORGANISM(S):

Escherichia coli ATCC 25922, Staphylococcus aureus ATCC 6538, Salmonella enterica subsp. Enterica ATCC 14028, Pseudomonas aeruginosa ATCC 9027, Streptococcus agalactiae ATCC 12386, Cronobacter sakazakii ATCC 29004, Listeria monocytogenes ATCC 7644

TEST RESULT(S):

As shown test result, comparing the level of viable organisms after sterilization to initial bioburden on the different treated feeding accessories, when using water volume of 100ml, this product can kill over 99.9999% microbe at the bioburden level of >6 log, either at 0 hour or 24 hour after sterilization.

Test microorganism: Escherichia coli ATCC 25922

Concentration of inoculated microorganism: 1.2×10^7 CFU/mL


The initial bioburden: 1.2×10^6 CFU/sample

Test time: The test was conducted at 0 hour post cycle completion.

Test Item	Test Group	Number of recovered microorganism (CFU/sample)	Efficacy(%)	Average Efficacy(%)	Log Reduction	Judge according to client's requirement: Log Reduction>6
Feeding bottle	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Feet	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Screw	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Cover	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Pump Valve	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			
Pump Body	1	<1	>99.9999	>99.9999	>6.0792	Pass
	2	<1	>99.9999			

SGS-CSTC Standards Technical Services (Shanghai) Co.,Ltd.

Page 2 of 9

Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request at www.sgs.com/portal/clients-english/conditions.aspx and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at File://www.sgs.com/portal/clients-english/conditions/Terms-Service-Documents.aspx. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained herein reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested.

3rd Building, No. 889, Yishan Road, Shanghai, China 200233

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Report No: ASH17-033560-03

Date: Dec 15 2017

Client name: Philips (China) Investment Co.,Ltd.
Client address: Philips Innovation Campus Shanghai No.1 Building, 10, Lane 888, Tian Lin Road Shanghai, P.R.C.
Sample name: Philips Avent 3-in-1 electric steam steriliser
Sample Batch No.: /
Product Date: /
Manufacturer: /
Model: SCF284

Above information and sample(s) was/were submitted and certified by the client, SGS quoted the information with no responsibility as to the accuracy, adequacy and/or completeness.

SGS Sample No.: ASH17-033560.001
Date of sample received: Aug 02 2017
Testing period: Aug 02 2017 ~ Aug 31 2017

TEST(S) REQUESTED:

Selected test(s) as requested by applicant
Validation Protocol of Electric Steam Steriliser Sterilization Effect

TEST METHOD(S):

Please refer to the next page(s)

TEST RESULT(S):

Please refer to the next page(s)

Chinese shall prevail in this report.


This Test Report supersedes the Test Report No. ASH17-033560-01 dated 19 Sep 2017 issued by SGS-CSTC Standards Technical Services (Shanghai) Co., Ltd. Original test report will be invalid from today.

Unless otherwise stated the results shown in this test report refer only to the sample(s) tested, and for clients internal use only, not to the society has the proof function. This document cannot be used for publicity, without prior written approval of the SGS.

SGS-CSTC Standards Technical Services (Shanghai) Co.,Ltd.

Page 1 of 9

Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com.cn/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com.cn/terms-and-conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, intent, location and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of a Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate or lie to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested.

3rd Building, No. 889, Yulin Road, Shanghai, China 200233

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www.sgs.com.cn

sgs.china@sgs.com

EUROPEAN DECLARATION OF CONFORMITY

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TUSSENDIEPEN 4, 9206 AD DRACHTEN, THE NETHERLANDS

(address / sede)

declare under our responsibility that the product(s): SCF284, SCF285

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to which this declaration relates is in conformity with the following harmonized standards:

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EC/1275/2008

2011/65/EU

And are produced under a quality scheme at least in conformity with ISO 9001 or CENELEC Permanent Documents

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Only for Medical Devices and R&TTE products:

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TUSSENDIEPEN 4, 9206 AD DRACHTEN, THE NETHERLANDS

(address / cím)

declare under our responsibility that the product(s) SCF284, SCF285

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